



Did you know?

As little as *seven minutes* of vigorous exercise can raise tendon and ligament temperatures to 113° F

How long do you exercise your horse?



Tendon Wraps



Evendura



Knee to Ankle
Wrap



Hock Wrap



Big Black Boot



Continuous Flow
Cold Water System



Ice Horse XC®

Immediate post-exercise cooling of the distal limbs of your horse can help prevent damaging, heat-related breakdown of tendon fiber.

Ice Horse products effectively cool targeted areas of your horse, snuffing out exercise-induced heat and inflammation *before* they become a significant threat to his soundness.